

POSITION GUIDE

ACHIEVABLE POSITIONS









Napping

Relaxing

Infinite Trendelenburg





2-POSITION TV RECLINE

Lifted • Seated • Reading

The 2-position TV recline is achievable on all Pride® Power Lift Recliners This is a comfortable position for performing many relaxing activities including reading, watching TV, working on your tablet or for conversation with family and friends.



3-POSITION **FULL RECLINE**

Lifted • Seated • Reading Napping

The 3-position full recline is available on most Pride® Power Lift Recliner models and puts you in a position ideal for napping or relaxing.



INFINITE-POSITION

Lifted • Seated • Reading Napping • Relaxing

The infinite-position recline is perfect for taking a nap. Infinite-position Pride® Power Lift Recliners enable the backrest and footrest to move independently from one another so that a variety of positions can be reached. The LC580i and LC525i can achieve true-infinite (two, three, infinite, zero gravity, Trendelenburg) positions while the LC558 can only achieve infinite (two, three, infinite) and lay flat positions.



TRENDELENBURG POSITION

The Trendelenburg position, available on the LC-525i and LC-580i, elevates one's feet higher than his or her head for increased comfort and relaxation.



ZERO GRAVITY

POSITION

Lifted • Seated • Reading • Napping • Relaxing Trendelenburg • Zero gravity

The zero gravity position, which can be reached on infinite-position Pride® Power Lift Recliner models, positions one's thighs and lower legs higher than his or her torso.